


LAB 1. IMAGE PRE-PROCESSING


Aim and data

2


- Understand the effects of each of the pre-processing steps seen in the class.
 - ▣ Bias field removal. Use MICO algorithm (matlab implementation).
 - ▣ Inter-patient normalisation
 - ▣ Anisotropic diffusion.
- DATA: 6 breast T2 MRI images from different patients (in png format)




pa1-16_t2.png




pa2-16_t2.png




pa3-16_t2.png



pa4-16_t2.png



pa5-16_t2.png



pa6-16_t2.png

1.1 Bias field removal (1h)

3

- Download the MICO algorithm for bias field removal.
- Investigate its parameters and how to run it with default images.
- Run the program for all data and generate the bias-free images.
- Deliver
 - ▣ Matlab code.
 - ▣ Image results before and after bias removal.

1.2 Patient Normalisation (1 h)

4

- Do the histogram of all patients
- Perform intensity normalisation with one of the methods seen in the class.
- Deliver
 - ▣ Matlab code.
 - ▣ Image results before and after normalisation.
 - ▣ Histogram before and after normalisation.
 - ▣ Bonus: compare without bias-field removed images.

1.3 Anisotropic diffusion (1 h)

5

- Evaluate the effects of each of the parameters of the anisotropic diffusion
 - ▣ Understand the file `anisoDiff2D.m`
 - ▣ See the effects of the parameters on the MRI and Ultrasound image.
 - ▣ Use `dicomread` and `dicominfo` to read the images.
 - ▣ Modify the code to perform isotropic smoothing (similar to Gaussian smoothing).
- Deliver
 - ▣ Matlab code.
 - ▣ Image results before and after smoothing.